

MAGENTA 2-WAY CAP

For Plant Tissue Culture and Test Tubes



FEATURES & BENEFITS:

- Drip Ring prevents condensation between the test tube and cap wall, minimizing contamination.
- Cap groove accommodates filter to eliminate contamination by purifying incoming air.
- Flexible ribs hold test tube snugly and compensate for outside diameter variations of test tubes.
- Visual transparency - natural or color tinted in blue, pink, yellow or white.
- Autoclavable and reusable, therefore reducing cost per use and unnecessary waste.
- Comfortable to handle.
- Lightweight.

COMMENTS

- 25 mm diameter x 31 mm height
- Works in conjunction with 25 mm Test Tubes
- Material: Polypropylene

ORDER INFORMATION

Part Number	Minimum Order Quantity	Pieces Per Box
2-Way Cap	1 Box	1,000 Pieces

HEADQUARTERS
15160 New Avenue
Lockport, Illinois 60441
TEL: (773) 777-5050

www.magentalc.com
customerfirst@magentalc.com

