


## MAGENTA 7-WAY TRAY

For Plant Tissue Culture and Test Tube Storage and Transport



### COMMENTS

- 11 7/8" L x 6" D x 1 5/8" T
- Material: Polypropylene

### FEATURES & BENEFITS:

- Minimized handling of tissue vessels and test tubes.
- Improves laboratory efficiency.
- Autoclavable and reusable, reducing costs and unnecessary waste.

### USES:

- **Way 1:** Carries vessels or culture tubes through continuous plant tissue culture operations including washing, filling, autoclaving, transporting and incubating without extra handling.
- **Way 2:** Holds (8) MAGENTA GA-7 or GA-7-3 Vessels.
- **Way 3:** Holds (36) 25mm culture tubes.
- **Way 4:** As a shelving module, resting on incubation rack channels instead of shelves.
- **Way 5:** Ventilates and eliminates hot spots between shelves, allows more closely spaced shelves on rack, reducing incubation costs.
- **Way 6:** Transmits light and allows more even lighting throughout the tray.
- **Way 7:** Lightweight and easily carried with one hand.

### ORDER INFORMATION

Part Number	Minimum Order Quantity	Pieces Per Box
7-Way Tray	1 Box	50 Pieces

**HEADQUARTERS**  
15160 New Avenue  
Lockport, Illinois 60441  
TEL: (773) 777-5050

[www.magentalc.com](http://www.magentalc.com)  
[customerfirst@magentalc.com](mailto:customerfirst@magentalc.com)

